


Setback and Clearance Notes

- 1 To facilitate service and filter changes, maintain a clearance of 28" from stern to nearest wall.
- 2 To maintain good air circulation for cooling, please maintain a minimum of 2" from the Starboard and Port sides to nearest wall.
- 3 Step access clearance 36" Minimum



Materials Data

Dry Weight: 463 lb
Water Weight: 1708 lbs
Wet Weight: 2231 lbs
Tub: 14GA 304 Stainless Steel - #4 Brushed
Staves: Incense Cedar
Woodframe: Incense Cedar
Steps: Incense Cedar
Wood Sealer: Exoshield Tung Oil

Electrical Data

Voltage: 120 volts single phase
Frequency: 60hz
Amps: 12.5 amps
GFCI Trip Current: 4-6mA

Mechanical Data

Compressor Power: 1/2 hp
Condensing Unit: AE 4450Y-AA1ADC Accumulator and Dryer Filter Attached
Refrigerant: R134a
Plunges Time per Day: 270 minutes
Ambient Operating Temperature: : 15°F - 120°F
Water Temperature: 0°F - 70°F

Water Treatment

Total Water Volume: 212 gal
Disinfection: O₃
Pump: DC80E-24100VR-1
Filter: 50sf Waterways, (Chlorine optional)
Ozone Injection Rate: 200 mg/hr (dual 100 mg/hr O₃ generator)
Flowrate: 9.8 GPM

DISCLAIMER:

THE MOROZKO IS DESIGNED TO SUPPORT A HEALTHY LIFESTYLE, NOT DIAGNOSE, CURE, OR PREVENT SPECIFIC DISEASES OR MEDICAL CONDITIONS. MOROZKO ICE BATHS ARE NOT MEDICAL DEVICES, AND HAVE NOT BEEN EVALUATED BY THE FDA. MOROZKO MAKES NO CLAIMS OR WARRANTIES WITH REGARD TO THE HEALTH EFFECTS OF COLD PLUNGE THERAPY OR USE OF MOROZKO PRODUCTS. SEEK MEDICAL ADVICE FROM YOUR PHYSICIAN BEFORE EMBARKING ON A PROGRAM OF DELIBERATE COLD EXPOSURE.